

LED OPERATION THEATER LIGHT

LED Operation Theater Light : 4 Clustor 4 Clustor 1,25,000 Lux



Model : CL4



Light Intensity	≤ 1,25,000/90,000 Lux
LED's Life	≤ 40,000 Hours
Color rendering Index	≤ 92
Temperature increase	≥ 2° C
Color Temperature	3500° – 6200°K
Input Supply	A C 220 – 230 V/50 Hz


Salient Features :

- Extremely environment friendly with least energy consumption
- Considerable longer life compared to Halogen bulb
- Combination of yellow & white LEDs for optimum comfort to Surgeon
- Shadow free uniform illumination inside cavities as well as open Operating area
- Control of Intensity by means of Digital Control Panel
- Easy adjustment and focusing of light by sterilizable handle
- Spring controlled self balancing arm with smooth radial, axial and angular movement.
- Suitable for air movement produced by laminar air flow.
- Most Comfortable for all major and multiple surgeries

LED Operation Theater Light : 3 Clustor 3 Clustor 90,000 Lux



Model : CL3



Model : M4



Model : M3



Combination of Dual and Triple Dome as well as Mobile – Stand Model in Any O T Light is available on special Order



ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgical@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice