

SUCTION MACHINE

Suction Machine - Mono




- Best Quality ¼ HP Single Phase Electric Motor
- Twin Bottle (2/1.5 Ltr) with overflow valve
- Self sealing rubber bush with nontoxic Rubber pipe
- Oil immersed Vacuum Pump
- Heavy duty castors for easy mobility

Suction Machine - Highvac



Suction Machine : M F A W 4L

- Works on main or by external 12 V D C Supply
- Built in rechargeable battery for emergency
- Collect up to 4 ltr. Fluids
- Polycarbonic jars with over flow valves
- Maintenance free pump due to oil free technology
- Vacuum up to 560 mm Hg with control
- Flow up to 22 LPM with bacterial filter
- Sp. Accessories : Foot Switch & 12 V Power Cable for use in vehicle
- Best for use in Ambulance



ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177



ACCUMAX

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice